



Fogg System Company, Inc.
15592 East Batavia Drive
Aurora, Colorado 80011
United States of America

Toll Free: 1 800 525-0292
Telephone: 303 344-1883
FAX: 303 344-1780

ECG WAVEFORM MEMORY MODULE INFORMATION

Part Number 7147-0041
DC Voltage Levels.

SHAPE Switch setting	AMPLITUDE ¹
AN	1000 μ V
AW	-1000 μ V
BN	Zero
BW	Zero

¹ *The amplitudes are nominal values at LA-RA for an AMPLITUDE setting of 1 mV using Fogg M310 & M311 Simulators, or equivalent.*